

PulmoAI: A Multi-Agent Clinical Decision Support System for Chest CT Analysis

Team: Med-Squad

- Ung Hoàng Long
- Trần Bảo Trân
- Đinh Trần Duy Trường
- Nguyễn Lê Thanh Minh

ABSTRACT

Problem

Lung cancer is the **#1 cause of cancer death** — yet diagnosis relies on scarce specialist expertise

Gap

Current AI diagnostic tools address only a single specialty, offering no integrated clinical view from a single scan.

Solution

PulmoAI: a multi-agent AI system that analyzes chest CT scans across 5 specialties in one pipeline

Method

MedGemma + LoRA · VoxTell auto-segmentation · LangGraph conversational agent

Deployment

Full-stack web app — accessible without specialist hardware

RELATED WORK

System	Specialty	System Capabilities
CheXNet	Radiology only	Radiology Text Only (No segmentation/dialogue)
Med-PaLM 2	General medical Q&A	Text Q&A Only (Not grounded in patient's scan)
LLaVA-Med	Medical VQA	Single-turn VQA (No clinical workflow)
nnU-Net / SAM-Med	Segmentation only	Pixel Segmentation Only (No report generation)
PulmoAI (Our)	Radiology + Cardiology + Oncology + Findings + Impression + Segmentation	Full End-to-End Pipeline

OBJECTIVES

- **Build a 5-agent AI pipeline for multi-specialty chest CT analysis**
- **Automate CT segmentation guided by the generated report**
- **Enable clinician-facing dialogue grounded in each patient's own results**
- **Deploy as a production web application — no specialist hardware required**

BUSINESS CONTEXT & DEPLOYMENT MODEL

Model	B2B — architecture-as-a-service for healthcare providers
How it works	Hospital provides their own clinical CT data → PulmoAI adapts the pipeline via fine-tuning
Why adaptable	Base model (MedGemma) is general-purpose — must be fine-tuned on hospital-specific data to reflect local patient population & terminology
This project	OpenM3Chest, CT-RATE used as simulated hospital dataset — MedGemma fine-tuned from scratch on it to demonstrate the adaptation pipeline

DATASET & PRE-PROCESSING & FINE-TUNING

DATASET OVERVIEW: OPENM3CHEST

Source	OpenM3Chest — derived from NLST (National Lung Screening Trial)
CT scan	Original DICOM (.dcm) → preprocessed to NumPy arrays in dataset
Labels & clinical data	JSON — demographics, smoking history, VQA questions, answers
Total tasks	15 labeled tasks across 3 clinical domains
Scope	Radiology · Cardiology · Oncology — no labels for Interpretation or Verification

OPENM3CHEST: TASK BREAKDOWN

Agent	Task	Output type
Radiology	Chest abnormality (×7)	Classification label
	Nodule presence	Yes / No
	Nodule attenuation	Solid / Ground-glass / Mixed
	Nodule location	Lobe (RUL, RLL, LUL...)
	Nodule margin	Spiculated / Smooth / Lobulated
	Nodule size	Measurement category
Cardiology	CVD diagnosis	Classification label
	CVD mortality risk	Low / High
Oncology	Lung cancer risk	No cancer / Cancer within follow-up

Radiology — Screening

Atelectasis

No

Is there any atelectasis, segmental or greater?

Pleural Effusion

Yes

Is there any pleural thickening or effusion?

Hilar / Mediastinal Mass

No

Is there a non-calcified mass or adenopathy $\geq 10\text{mm}$ in the hilar/mediastinal region?

Chest Wall Abnormality

No

Is there any chest wall abnormality (bone destruction or metastasis)?

Consolidation

No

Is there any consolidation?

Emphysema

No

Is there any emphysema?

Fibrosis / Honeycombing

No

Is there any fibrosis, honeycombing, reticular/reticulonodular opacities, or scarring?

Nodule Present

Yes

Is there any lung nodule?

Radiology — Nodule Detail

Size

4~6mm

What is the nodule size?

Margin

Smooth

What is the nodule margin type?

Location

Right Lower Lobe

In which lobe is the nodule located?

Attenuation

Solid

What is the predominant attenuation of the nodule?

Cardiology

CVD Diagnosis

Yes

Can you identify any notable abnormalities in the cardiovascular system?

CVD Mortality Risk

Low risk

Predict the risk of cardiovascular disease mortality.


Oncology

Lung Cancer Risk

No cancer within follow-up

Given the current indicators, what's the six-year outlook for developing lung cancer?

Data Preprocessing & Parallel Multi-Adapter LoRA Fine-Tuning



DATASET OVERVIEW: CT-RATE

Source	CT-RATE: Open source dataset for 3D CT-SCAN
CT scan	Original DICOM (.dcm) → preprocessed to NumPy arrays in dataset
Total tasks	2 Tasks: Findings and Impressions


NOTES:

- **Findings:** What the radiologist **sees** (observations, measurements, abnormalities on the scan).
- **Impressions:** What those findings **mean** (summary, final diagnosis, or next steps).

DATASET OVERVIEW: CT-RATE

Why CT-RATE over OPENM3CHEST?

- **Limitation of OPENM3CHEST:** Restricted to simple binary classification (Yes/No labels).
- **The CT-RATE Advantage:** Provides rich, narrative data structured into clinical Findings and Impressions.
- **Clinical Impact:** Moves the AI from a basic classifier to a true diagnostic assistant, generating detailed reports that actively support radiologists in complex decision-making.



DATASET SPLIT FOR AGENT FINE-TUNING


	OPENM3CHEST	CT-RATE
Train	Oncology: 500 Radiology: 500 Cardiology: 500	Findings and Impressions: 500
Test	Oncology: 200 Radiology: 200 Cardiology: 200	Findings and Impressions: 200

 **Small subset used due to GPU resource constraints (A100 80GB required)**

METHOD

Multi-Agent Clinical Report Generation Pipeline

MULTI-AGENT CLINICAL REPORT GENERATION PIPELINE



VERIFICATION AGENT

What it does

- Independently re-reads the raw report from all 4 agents
- Can disagree with specific findings
- Surfaces unreported findings missed by specialist agents
- Adds clinical caveats and differential considerations

How it works

- Uses base MedGemma — no LoRA adapter
- Intentionally unbiased: not fine-tuned toward any specialty

Output

- Verified report → passed to Highlight agent (DeepSeek) for color annotation

FINDINGS AND IMPRESSION

Chest

I agree with all points regarding the **absence of atelectasis**, **pleural thickening/effusion**, **consolidations**, parenchymal lesions like **masses or nodules >1cm**, **bone destruction**, etc., as well as the lack of obvious signs suggestive of **severe pulmonary fibrosis** such as extensive honeycombing or widespread reticular patterns beyond mild background interstitial markings potentially related to chronic inflammation from her known respiratory conditions. **Emphysematous changes appear prominent bilaterally** involving primarily the apices but also extending into more basal regions. Some **patchy groundglass opacity** might represent superimposed inflammatory processes rather than definite edema. **Small calcifications along bronchovascular bundles** may reflect benign granulomas. I **disagree with the finding "Is there any consolidation?"** → Yes. While some **subtle hazy increased density** exists particularly towards the right middle lobe base possibly representing **early infectious process** versus postinflammatory changes, definitive lobar consolidation isn't clearly visualized across multiple phases if available. However, small foci cannot entirely rule out **very early pneumonia** especially considering patient symptoms. If clinical suspicion remains high despite negative imaging, further evaluation including lateral decubitus views or short interval repeat exam would help clarify. Also, consider correlation with laboratory data.

Cardiovascular

I agree with the assessment showing **unremarkable cardiac silhouette** without **pericardial fluid accumulation** nor **aortic dissection features**. **Normal caliber vessels** throughout.

Oncology Risk

Based solely on these static images presented here, while she has several factors increasing potential future malignancy risk due to **heavy smoking history** and **possible underlying COPD** contributing to impaired clearance mechanisms, currently nothing specific suggests imminent development of primary lung cancer compared to general population risks unless additional information becomes available suggesting otherwise. This interpretation should always account for limitations inherent in assessing long term prognosis purely through single snapshot scans.

Imaging Report Accuracy

I largely agree with the radiology summary section describing the overall picture including **emphysema pattern**, **bibasilar linear scar tissue** likely reflecting old pneumonic infiltrates leading to atelectasis ("sequelae density increase"), and the **solitary tiny (<6mm) solid nodule noted in the left lower lobe periphery**. It accurately reflects common incidental findings often seen in patients with advanced obstructive airways diseases who smoke heavily.


Unreported Findings

Besides the aforementioned minor considerations about possibility of resolving infection vs postinflammation causing slight haziness near RML base, another observation includes **scattered punctate micronodules distributed diffusely throughout both lungs measuring up to ~3mm** maximum dimension. These can sometimes be associated with various etiologies ranging from simple intrapulmonary hemorrhage secondary to vasculitis affecting smaller vessels, metastatic seeding although less likely given their appearance, or even reactive lymphoid hyperplasia. Their significance depends greatly upon context and prior medical history. Another noteworthy feature is **marked tortuosity and atherosclerotic plaque burden evident in descending thoracic aorta** and proximal abdominal segments shown partially imaged. **Mild degenerative discogenic spurring** is apparent lumbar spine levels.

FUNCTION: CHATBOT

FUNCTION: CHATBOT

CHAT MEMORY



PostgreSQL




CLINICAL DATA OF
PATIENT



REPORT OF CT-SCAN



MCP



- **Drug Interactions Tool:** Identifies risk factors when mixing medications.
- **Warnings & Cautions Tool:** Extracts safety alerts, contraindications, and clinical precautions.

FUNCTION: CHATBOT

Drug Interaction

Ask AI

The doctor prescribed Symbicort (or Spiriva) for this. Is this medication safe and optimal for his respiratory profile?"

Now I have comprehensive data from the drug safety profiles and interaction databases. Here is my detailed evidence-based analysis.

Clinical Analysis: Safety & Optimality of Symbicort vs. Spiriva for This COPD Patient

Patient Profile Summary

Risk Factor	Details
Age	63 years
Sex	Female
BMI	18.9 (Underweight)
COPD severity	Bilateral emphysema, chronic bronchitis since childhood
Smoking	90 pack-years (former), lives with smoker
Prior pneumonia	Yes (age 36)
Cancer Hx	Breast cancer (age 44)

Ask about this scan...

Summarization

Can you summarize the report of the patient for me please?

Certainly! Here is a clear and organized summary of your medical report:

Patient Summary

Patient: 63-year-old White female

Height: 61 inches (5'1") | Weight: 100 lbs (BMI 18.9)

Smoking: Former heavy smoker (90 pack-years), quit at age 51, currently lives with a smoker.

Key Medical History:

Chronic bronchitis (since age 4)

COPD (since age 47)

Pneumonia (since age 36)

Breast cancer (age 44 – in remission)

Family history: Mother had lung cancer

Ask about this scan...

FUNCTION: SEGMENTATION


FUNCTION: SEGMENTATION

System Overview:

- The VoxTell Web system allows doctors to view CT scans on a web interface and automatically generate an initial segmentation mask based on the patient report.


Component	Role
Frontend Vite	Displays CT volume and segmentation masks
FastAPI Backend	Handles API requests and coordinates the workflow
Modal GPU	Runs VoxTell segmentation with GPU acceleration
S3 + Postgres	Stores DICOM data and patient analysis reports

Overall Architecture



i Local deployment: Backend proxies requests to Modal via `VOXTELL_MODAL_BASE_URL`

FUNCTION: SEGMENTATION



FUNCTION: SEGMENTATION

**VoxTell CT Viewer**

AUTO DICOM FROM S3 · NIIVUE VOLUME VIEWER

done




Loading CT volume from S3...
pid=122833 · series_uid=1.2.840.113654.2.55.307229896835702397101325832166669748691

AUTO SEGMENTATION FROM REPORT

Hệ thống sẽ tự đọc report, dùng LLM chọn target quan trọng, thử nhiều prompt và chọn mask tốt nhất.


MANUAL TEXT PROMPT

VD: left ventricle, tumor, lung nodule...



VoxTell CT Viewer

Waiting for backend volume conversion





No CT volume loaded
Backend sẽ tự tải dicom-raw/{pid}/{series_uid}.zip từ S3 và convert sang NIfTI tạm.

SEGMENTATIONS


0

Chưa có mask segmentation. Auto mask sẽ xuất hiện trước; bác sĩ vẫn có thể nhập prompt thủ công để thêm mask mới.

FUNCTION: SEGMENTATION



Auto Segmentation Workflow



Example

- Top-3 findings
 - lung nodule
 - pulmonary nodule
 - right lower lobe lung nodule

Selected prompt: **lung nodule**

Validation: **nonzero voxels > 0**




Auto segmentation provides an initial suggestion. Doctors can still run **manual prompts**.

FUNCTION: SEGMENTATION



FUNCTION: SEGMENTATION



FUNCTION: SEGMENTATION

**VoxTell CT Viewer**

AUTO DICOM FROM S3 · NIIVUE VOLUME VIEWER

done

**CT volume ready**

pid=122833 · series_uid=1.2.840.113654.2.55.307229896835702397101325832166669748691

AUTO SEGMENTATION FROM REPORT

Auto segmentation đã sẵn sàng với prompt tốt nhất.

TOP-3 SELECTED BY LLM

1. Emphysematous changes with peripheral reticulation in both lungs

2. 3 mm nonspecific nodule in left lower lobe superior segment laterobasal area

3. Sequelae fibrotic densities compatible with atelectatic change in the apex of both lungs

Selected VoxTell prompt: lung nodule

PROMPT CANDIDATES

☒ lung nodule nonzero=1352

MANUAL TEXT PROMPT

heart

Processing...

VoxTell CT Viewer

See more clearlyReset

Current view: Multi


S

L

P

A

L



WL

-959..437

Axial

Coronal

Sagittal

Multi

SEGMENTATIONS

2

AUTO

lung nodule

Top-3: Emphysematous changes with peripheral reticulation in both...

MANUAL

heart

FUNCTION: SPEECH TO TEXT

FUNCTION: SPEECH TO TEXT

Purpose:

- Allow doctors to record their clinical comments or questions directly via a Microphone button in the browser.
- Automatically transcribe speech to specialized English medical text in real-time and insert the result into the input field for medical Chatbot queries.
- Help doctors minimize manual typing time during complex case consultations

FUNCTION: SPEECH TO TEXT

Model Specification:

- core model: OpenAI Whisper-Base model (73 million parameters) fine-tuned on a specialized medical terminology dataset focusing on Medicine and Surgery
- Inference Library: faster-whisper
- Resources Allocated: Runs entirely on CPU (cpu=2.0 - 2 vCPUs, memory=2048 - 2GB RAM).

RESULTS

RESULTS

Evaluation Metrics for Clinical Report Generation

Lexical Overlap Metrics (BLEU & ROUGE):

- **BLEU (1-4):** Measures n-gram precision. It checks how many words/phrases in the generated report exactly match the ground truth.
- **ROUGE (1, L):** Measures recall. It ensures the model doesn't miss key clinical information from the reference.

Semantic Similarity (BERT-Score):

- **The "Gold Standard":** Unlike BLEU/ROUGE, BERT-Score uses contextual embeddings to measure meaning.
- **Why it matters:** In medicine, "Heart is enlarged" and "Cardiomegaly is present" mean the same thing. BERT-Score recognizes this similarity even if the words are different.

Classification Accuracy:

- Used for specific diagnostic tasks (Oncology/Cardiology/Radiology) to measure the percentage of correct clinical predictions.

RESULTS

Findings and Impressions Agent

Model	Task	BLEU	BLEU-1	BLEU-4	ROUGE-1	ROUGE-2	ROUGE-L	BERT-SCORE
Finetuned	Findings	0.0652	0.1121	0.0447	0.3815	0.2568	0.3015	0.9398
	Impressions	0.0307	0.2938	0.0397	0.1702	0.1073	0.1503	0.9147
Base	Findings	0.0033	0.0337	0.0003	0.1595	0.0260	0.1037	0.9126
	Impressions	0.0000	0.2514	0.0000	0.0833	0.0159	0.0740	0.9151

- **Vocabulary & Structure (BLEU/ROUGE):** The Base model struggles heavily with medical reporting standards (BLEU-4 is **nearly 0.0000**). Fine-tuning successfully forces the model to adopt the correct terminology and narrative flow, multiplying ROUGE-1 and ROUGE-L scores significantly.
- **Clinical Meaning (BERT-Score):** Interestingly, the Base model achieves a high BERT-Score (~0.915). This indicates the base LLM correctly identifies the underlying pathologies, but simply lacks the domain-specific vocabulary to express them.
- **The Fine-Tuning Impact:** LoRA fine-tuning effectively acts as a "formatting and vocabulary bridge," teaching the AI how to write like a real radiologist without losing its core semantic intelligence.

RESULTS

Oncology Agent

Task / F1-Score	BASE (No Lora)	NO EXPLAIN (Finetuned)	WITH EXPLAIN (Finetuned + CoT)
Lung cancer risk	0.52	0.5	0.46

Slight Regression

- Base accuracy experienced a minor drop (from **0.52** to **0.46** using CoT).
- **Conclusion:** Suggests potential localized overfitting, imbalanced training data, or noise introduced by the Chain-of-Thought prompting in this specific domain.

Cardiology Agent

Task / F1-Score	BASE (No Lora)	NO EXPLAIN (Finetuned)	WITH EXPLAIN (Finetuned + CoT)
CVD diagnosis	0.21	0.5	0.55
CVD mortality	0.1	0.94	0.78

Massive Improvement

- Successfully overcame baseline limitations.
- **Accuracy** for CVD mortality surged dramatically (from **0.10** to **0.94**).
- **Conclusion:** The model effectively learned complex cardiovascular clinical patterns.

RESULTS

Radiology Agent

Task / F1-Score	BASE (No Lora)	NO EXPLAIN (Finetuned)	WITH EXPLAIN (Finetuned + CoT)
Chest abn 54	0.52	1	1
Chest abn 55	0.58	0.21	1
Chest abn 56	0.59	0.86	0.97
Chest abn 57	0.61	1	1
Chest abn 58	0.58	0.97	0.97
Chest abn 59	0.47	0.59	0.51
Chest abn 61	0.39	0.7	0.65
Nodule attenuation	0.06	0.27	0.18
Nodule location	0.06	0.27	0.18
Nodule margin	0.26	0.73	0.73
Nodule presence	0.29	0.58	0.51
Nodule size	0.04	0.38	0.38

RESULTS

Radiology Key Findings

- **Essential Specialization:** The Base model failed significantly in Nodule tasks (~0.04 - 0.06 accuracy). Fine-tuning was mandatory to achieve diagnostic capability.
- **The CoT Breakthrough:** Chain-of-Thought (With Explain) achieved Perfect Accuracy (1.0) across multiple Chest Abnormalities (54, 55, 57).
- **Reasoning Recovery:** Notably in Chest abn 55, the direct fine-tuned model dropped to 0.21, but CoT recovered it to 1.0, proving that explicit reasoning prevents logical collapses.
- **Pattern vs. Logic:** No-Explain performed better on visual-spatial tasks (Nodules), while CoT excelled in complex diagnostic logic.

The Power of Fine-Tuning

- **Fine-Tuning > Base Model:** Across all domains (Oncology, Cardiology, Radiology), fine-tuning transformed a non-functional base model into a specialized medical expert.
- **Clinical Viability:** While Base models struggle with zero-shot medical reasoning, LoRA fine-tuning provides the domain-specific knowledge necessary for real-world clinical use.
- **Final Verdict:** Despite minor variances in Oncology, the massive accuracy gains (up to 9x improvement) in specialized tasks prove that domain-specific fine-tuning is the only path to reliable AI diagnostics.

DISCUSSION

CLINICAL VALUE OF MULTI-AGENT DISCREPANCIES

How Agent Conflicts Empower Clinicians

1. Eliminating AI Confirmation Bias

- Single-model AI often forces doctors into "blind trust."
- In our system, discrepancies between the Base Agents and the Verification Agent serve as an automated clinical red flag, preventing physicians from overlooking critical omissions.

2. Active Human-in-the-Loop Triage

- When agents disagree, the system does not make an arbitrary choice. Instead, it flags the case as a High-Risk Mismatch.
- This mismatch instantly alerts the doctor to actively intervene and utilize Feature 2 (Interactive Segmentation) or Feature 3 (RAG Chat) to audit the exact medical pixels.

3. Multi-Perspective Diagnostic Safety


- By exposing the internal conflict (e.g., text vs. image, focal vs. diffuse), the framework provides doctors with a comprehensive, multi-angle view of the patient's condition, converting technical boundaries into clinical safety nets.

CLINICAL EVALUATION - CASE 1: LUNG NODULE

COMPONENT	RADIOLOGY AGENT	VERIFICATION AGENT	GROUND TRUTH
Lung Nodule Detail	• Location: Right Lower Lobe • Size: 4~6mm • Margin: Smooth	• Correction: Relocated to the Left Lung • Refinement: Solitary nodule is < 6mm (~3mm)	• Location: Left Upper Lobe • Size: 6~8mm • Margin: Poorly defined

The Conflict: Base agents flag a smooth nodule in the Right Lower Lobe, while the Verification agent identifies it in the Left Lung. Meanwhile, Ground Truth marks a dangerous Poorly-defined nodule in the Left Upper Lobe.

“Segment the poorly defined solid nodule in the left upper lobe.”



CONCLUSION

CONCLUSION

- This project successfully develops a web-based medical AI system for supporting chest CT analysis.
- The system integrates DICOM preprocessing, clinical data processing, LoRA fine-tuned agents, report generation, verification, and lesion segmentation.
- Experimental results show that LoRA fine-tuning significantly improves model performance, especially for the Radiology Agent, Cardiology Agent, and Findings and Impressions Agent.
- Fine-tuning not only improves prediction performance but also helps the model generate reports with more appropriate medical terminology and structure.

REFERENCES

- [1] I. E. Hamamci et al., "Generalist foundation models from a multimodal dataset for 3D computed tomography," *Nature Biomedical Engineering*, pp. 1–19, 2026.
- [2] J. Hofmanninger, F. Prayer, J. Pan, S. Röhrich, H. Prosch, and G. Langs, "Automatic lung segmentation in routine imaging is primarily a data diversity problem, not a methodology problem," *European Radiology Experimental*, vol. 4, no. 1, p. 50, 2020.
- [3] A. SELLERGRÉN et al., "MedGemma 1.5 technical report," *arXiv preprint arXiv:2604.05081*, 2026.
- [4] C. Niu, Q. Lyu, C. D. Carothers, P. Kaviani, J. Tan, P. Yan, M. K. Kalra, C. T. Whitlow, and G. Wang, "Medical multimodal multitask foundation model for lung cancer screening," *Nature Communications*, vol. 16, no. 1, p. 1523, 2025.
- [5] J. Ma, Y. He, F. Li, L. Han, C. You, and B. Wang, "Segment anything in medical images," *Nature Communications*, vol. 15, no. 1, p. 654, 2024.
- [6] M. Rokuss et al., "VoxTell: Free-text promptable universal 3D medical image segmentation," in *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 2026, pp. 37538–37557.

TECHSTACK




CLOUD & SERVERLESS



WEB APPLICATION

WEB CLIENT




BACKEND



AI



unsloth



LangChain

DATABASE




PostgreSQL



Amazon S3

EXTERNAL SERVICE



deepseek




HUGGING FACE




GitHub

DATABASE DESIGN

1. POSTGRESQL REAL DATABASE SCHEMA



2. AWS S3 OBJECT STORAGE





THANK YOU

FOR LISTENING

End Slide